

Ambu® Redi-ACE Adult & Redi-ACE Mini



Key Benefits

- **Easy to use**
Ambu Redi-ACE has an automatic flip up chin and two locking buttons. It can be placed on the patient using only 3 steps; sizing – locking – placement
- **14 Sizes in 2 collars**
With 8 sizes of Ambu Redi-ACE Adult and 6 sizes of Ambu Redi-ACE Mini the majority of patients can be treated
- **Practical characteristics**
Ambu Redi-ACE can be stored completely flat. Additionally the Ambu Redi-ACE has tube hooks to secure the oxygen tube
- **CT, MR and X-ray compatible***

Ambu® Redi-ACE Adult & Redi-ACE Mini

The Ambu Redi-ACE is a unique adjustable collar. Combined with the Redi-ACE Mini these two collars offer 14 sizes ensuring that patients can have the most effective and accurate fit possible.

The collar is packaged completely flat and has a unique automatic flip up chin.

Ambu Redi-ACE is a disposable adjustable collar designed to assist with the maintenance of neutral alignment, prevention of lateral sway and anterior-posterior flexion and extension of the spine during transport.

Two locking buttons hold the collar in place when adjusted.

* Distortion or shadows may appear on the image depending on the settings of the imaging equipment

Specifications

Ambu Redi-ACE	
Adjustable sizes	Neckless, Short, Regular and Tall
Number of settings	8 adjustments
Dimensions	554 x 215 x 14 mm (21.8" x 8.5" x 0.6")
Weight	200 g (7 ounces)
Color	Blue/light grey
Translucency	X-ray, CT and MRI compatible
Storage temperature	-30° to +70° C
Operation temperature	-10° to +40° C
Ambu Redi-ACE Mini	
Adjustable sizes	Small pediatric and Small Adult
Number of settings	6 adjustments
Dimensions	447 x 188 x 14 mm (17.6" x 7.4" x 0.6")
Weight	100 g (3.5 ounces)
Color	Yellow/light grey
Translucency	X-ray, CT and MRI compatible
Storage temperature	-30° to +70° C
Operation temperature	-10° to +40° C



Ambu Redi-ACE Adult



Ambu Redi-ACE Mini

Material

Product	Material
Ambu Redi-ACE Adult	
Blue/light grey plastic	Polypropylene
Foam	IXPE
Ambu Redi-ACE Mini	
Yellow/light grey plastic	Polypropylene
Foam	IXPE



Flat storage

Ordering specifications

Item no.	Description
472 001 000	Redi-ACE Adult Extrication Collar, 30/Case
472 002 000	Redi-ACE Mini Extrication Collar, 30/Case



Locking buttons



Tube hooks

Ambu A/S

Baltorpbakken 13
 DK - 2750 Ballerup
 Denmark
 Tel. +45 7225 2000
 ambu.com

